

## Glass DC Nebulizers

### SeaSpray

400µL/min uptake      A13-1-USS04X

### MicroMist

400µL/min uptake      A13-1-UM04X

200µL/min uptake      A13-1-UM02X

100µL/min uptake      A13-1-UM01Q

50µL/min uptake      A13-1-UM005Q



### Conikal

1mL/min uptake      A13-1-UC1X

## Quartz DC Nebulizers

### SeaSpray

1mL/min uptake      A13-1-QSS1



## DC Nebulizers for HF

### OpalMist

400µL/min uptake      A13-1-PFA04X

200µL/min uptake      A13-1-PFA02X

100µL/min uptake      A13-1-PFA01Q

50µL/min uptake      A13-1-PFA005Q



## Peltier Controlled Cyclonic Spray Chambers

Glass Spray Chamber      KT-1116

Quartz Spray Chamber      KT-1116Q

PFA Spray Chamber      KT-1116P



## Spray Chambers

### Single-Cell Spray Chamber Kit

High Sensitivity Single-Cell Sample Introduction System      KT-1155



### HydraMist Simultaneous Cold Vapor/Pneumatic Nebulization Spray Chamber

Hydride Kit      KT-1168



## Spray Chambers

### \*Twister

20mL Cyclonic with Aux.      20-809-2574HE

### \*Cinnabar

20mL Cyclonic with Aux.      20-809-2575HE



## Spray Chambers for HF

### \*Tracey PFA44

44ml Cyclonic with Aux.      20-809-2578

Spray Chamber Adaptor      KT-1104

\*Requires adaptor KT-1104.



## Torches


### Quartz Torches

with 1.0mm Injector      30-807-0558

with 1.5mm Injector      30-807-0557

with 2.0mm Injector      30-807-3608

with 2.5mm Injector      30-807-0556



### Torch Options

Quartz Torch Bonnet      31-808-2577

Quartz Connector with Dilution Port      31-808-3217

Quartz Connector with 45deg Dilution Port      31-808-3280

Straight Quartz Connector      31-808-3288

Clamp (PKT 2)      70-803-0846



### Semi-Demountable Torch

with 2.4mm Quartz Injector      30-808-3853



### Torch Options

2.4mm Tapered Quartz Injector      31-808-4168

2.4mm Tapered Alumina Injector      31-808-3855

2.4mm Sapphire Capillary Injector      31-808-4159



**RF Coils**

Copper	70-900-7800C	
Silver	70-900-7800S	
Gold	70-900-7800G	
<b>RF Coil Installation Kit</b>	70-900-7500T	

**Cones**

**Sampler Cones**

Nickel Sampler	AT7701-Ni	
Nickel Sampler with Nickel plated base	AT7701-Ni/Ni	
Aluminium Sampler	AT7701-Al	
Platinum Sampler	AT7706-Pt	
Platinum Sampler (18 mm Insert)	AT7706A-Pt	
Platinum Sampler Boron Free	AT7706-Pt-BF	

**Cones**

(continued)

**Sampler Cones (continued)**


Platinum Sampler with Nickel plated base	AT7706-Pt/NiP	
Platinum Sampler with solid Nickel base	AT7706-Pt/Ni	
Platinum Sampler (15mm insert) with Nickel base	AT7706B-Pt/Ni	
Platinum Sampler with Nickel plated base (Extended Life)	AT7706L-Pt/NiP	

**Skimmer Cones**

Nickel Skimmer for x-lens	AT7702X-Ni	
Aluminium Skimmer for x-lens	AT7702X-Al	
Platinum Skimmer for x-lens	AT7708X-Pt	

Platinum Skimmer with Nickel base for x-lens	AT7708X-Pt/Ni
Platinum Skimmer with Boron Free for x-lens	AT7708X-Pt-BF

**Cone Options**

Skimmer Cone Guard	70-803-1008	
Graphite Gasket (PKT 3)	AT7703	
Long Life Shield Plate	AT5004	

**Common Accessories**

**Elegra Argon Humidifier**

Elegra Dual Two-channel Argon Humidifier 70-803-1273



Elegra Stand 70-803-1581



**Guardian In-Line Filters**

Guardian DC in-line non-return gas filter 70-803-1942


Guardian in-line sample filter 70-803-1108



**Trident Internal Standard Kits**

Trident CT Standard Kit 60-703-1179

Standard Kit for HF 60-808-1150



**Autosampler Probes**

(See page 14 for complete list of probes)

**For Agilent® I-AS**

PTFE Encapsulated Carbon Fibre Probe 0.25mm ID 60-703-1010

PTFE Encapsulated Carbon Fibre Probe 0.3mm ID 60-703-1009

**For Cetac ASX-200/500/800 Series**

Polyimide sheathed Autosampler Probe 0.5mm ID 60-808-1186L

PTFE Encapsulated Carbon Fibre Probe 0.5mm ID 70-803-0784

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-1443

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0793



**For SPS 3/SPS 4**

PTFE Encapsulated Carbon Fibre Probe 0.75mm ID 70-803-0908

PTFE Encapsulated Carbon Fibre Probe 1.0mm ID 70-803-0853