


## Sampler cones for 7700/8800/7800/7900/8900 ICP-MS

Description	Part Number	Comments
Nickel sampler cone	AT7701-Ni	Standard with x-lens
Nickel sampler cone with Nickel plated base	AT7701-Ni/Ni	For samples with > 0.5% HCl
Aluminium sampler	AT7701-Al	
Platinum sampler cone	AT7706-Pt	Standard with s-lens
Platinum sampler cone (18 mm insert)	AT7706A-Pt	For use with high viscosity, high boiling point acids
Platinum sampler cone Boron Free	AT7706-Pt-BF	
Platinum sampler cone with Nickel plated base	AT7706-Pt/NiP	
Platinum sampler cone with Nickel base	AT7706-Pt/Ni	
Platinum sampler cone (15 mm insert) with Nickel base	AT7706B-Pt/Ni	



AT7701-Ni



AT7701-Ni/Ni

## Skimmer cones for 7700/7800/8800 ICP-MS

Description	Part Number	Comments
Nickel skimmer cone for x-lens	AT7702X-Ni	Standard with x-lens
Aluminium skimmer cone for x-lens	AT7702X-Al	
Platinum skimmer cone for x-lens	AT7708X-Pt	Best LOD in high matrix
Platinum skimmer cone with nickel base for x-lens	AT7708X-Pt/Ni	Recommended for organic analysis
Platinum skimmer cone Boron Free for x-lens	AT7708X-Pt-BF	
Nickel skimmer cone for s-lens	AT7702S-Ni	
Platinum skimmer cone for s-lens	AT7708S-Pt	Standard with s-lens
Platinum skimmer cone with nickel base for s-lens	AT7708S-Pt/Ni	



AT7708X-Pt



AT7708X-Pt/Ni

## Skimmer cones for 7900 ICP-MS

Description	Part Number	Comments
Nickel skimmer cone for x-lens	AT7902X-Ni	Standard with x-lens
Aluminium skimmer cone for x-lens	AT7902X-Al	
Platinum skimmer cone for x-lens	AT7908X-Pt	Best LOD in high matrix
Platinum skimmer cone with nickel base for x-lens	AT7908X-Pt/Ni	Recommended for organic analysis
Nickel skimmer cone for s-lens	AT7702S-Ni	
Copper skimmer cone	AT7702S-Cu	
Platinum skimmer cone for s-lens	AT7708S-Pt	Standard with s-lens
Platinum skimmer cone with nickel base for s-lens	AT7708S-Pt/Ni	



AT7902X-Ni



AT7908X-Pt

## Skimmer cones for 8900 ICP-MS

Description	Part Number	Comments
Nickel skimmer cone for x-lens	AT7902X-Ni	Standard with x-lens
Aluminium skimmer cone for x-lens	AT7902X-Al	
Platinum skimmer cone for x-lens	AT7908X-Pt	Best LOD in high matrix
Platinum skimmer cone with nickel base for x-lens	AT7908X-Pt/Ni	Recommended for organic analysis
Nickel skimmer cone for s-lens	AT8902S-Ni	
Copper skimmer cone for s-lens	AT8902S-Cu	
Platinum skimmer cone for s-lens	AT8908S-Pt/Cu	Standard with s-lens
Platinum skimmer cone with nickel base for s-lens	AT8908S-Pt/Ni	
Platinum skimmer cone with nickel base for m-lens	AT8908M-Pt/Ni	For use with m-lens



AT7908X-Pt/Ni



AT7902X-Ni

## Cone Options


Description	Part Number	Comments
Skimmer Cone Guard	70-803-1004	For 7700s/7900
Skimmer Cone Guard	70-803-1008	For 7700x/7800/8800/8900
Retaining Ring	AT7704	For 7700 Sampler Cone
Graphite Gasket (PKT 3)	AT7703	For 7700/7900/8800/8900 Sampler Cone
Long Life Shield Plate	AT5004	For 7700/7800/7900/8800/8900



70-803-1004



AT7704



AT5004