

ICP-OES/ICP-MS SUPPLIES & ACCESSORIES

Issue 9

Product List for Shimadzu ICP-OES/MS: 1000, 5000, 7000, 8000, & 9000



GLASS EXPANSION
Quality By Design

INTERNATIONAL

15 Batman Street
West Melbourne
Vic 3003, Australia

Telephone: +61 3 9320 1111
Facsimile: +61 3 9320 1112
Email: enquiries@geicp.com

AMERICAS

4 Barlows Landing Road, Unit 2A
Pocasset, MA 02559, USA

Toll Free Phone: 800 208 0097
Telephone: 508 563 1800
Facsimile: 508 563 1802
Email: geusa@geicp.com

Glass Nebulizers

Conikal
1mL/min uptake [ARG-1-JC1](#)

SeaSpray
2mL/min uptake [ARG-1-USS2](#)

Slurry
4mL/min uptake [ARG-1-US6](#)

MicroMist
400µL/min uptake [ARG-1-UM04](#)
200µL/min uptake [ARG-1-UM02](#)
100µL/min uptake [ARG-1-UM01](#)
50µL/min uptake [ARG-1-UM005](#)



Nebulizers suitable for use with HF

OpalMist
2mL/min uptake [ARG-1-PFA2](#)
400µL/min uptake [ARG-1-PFA04](#)
200µL/min uptake [ARG-1-PFA02](#)
100µL/min uptake [ARG-1-PFA01](#)
50µL/min uptake [ARG-1-PFA005](#)



DuraMist
1mL/min uptake [ARG-1-DM1](#)
400µL/min uptake [ARG-1-DM04](#)



VeeSpray - Ceramic
Standard [ARG-1-CV6](#)



IsoMist

IsoMist Programmable Temperature Spray Chamber for Shimadzu

With glass spray chamber [KT-1067](#)
With quartz spray chamber [KT-1067Q](#)
With PFA spray chamber [KT-1067P](#)



Spray Chambers


Tracey
50mL Cyclonic [20-808-8881HE](#)
With B14 gravity drain.



Twister
50mL Cyclonic [20-809-9119HE](#)
With B14 gravity drain.



Twister
50mL Cyclonic [20-809-0222HE](#)



Cinnabar
20mL Cyclonic [20-809-0164HE](#)



Drain Device
[21-808-8214](#)



Glass Transfer Tube, Spray Chamber to ICPE-9000

[21-809-2846](#)



Tracey TFE
50mL Cyclonic [20-809-2506](#)
Suitable for use with HF.



Spray Chambers (continued)

Tracey TFE

50mL Cyclonic [20-809-2508](#)
 With B14 gravity drain.
 Suitable for use with HF.



HF Drain Device

[21-809-2540](#)



Mounting bracket for HF Spray Chamber

[21-809-2839](#)



Torches (continued)

Quartz Torches

Quartz Torch Standard
 1.5mm Injector (Suitable for
 ICPS-7500/8100/9000) [30-808-1179](#)



Quartz Torch for Organics,
 1.0mm Injector (Suitable for
 ICPS-7500/8100/9000) [30-807-0519](#)

Quartz Torch for High TDS,
 1.8mm Injector (Suitable for
 ICPS-7500/8100/9000) [30-807-0518](#)

Quartz Torch 1.2mm Injector
 (Suitable for ICPM-8500) [30-808-9986](#)



Quartz Torch 1.5mm Injector
 (Suitable for ICPM-8500) [30-808-1090](#)

Quartz Torch L Shape
 (Suitable for ICPS-7000) [30-808-1091](#)



Quartz Mini Torch
 (Suitable for ICPE-9000) [30-807-0534](#)



Torches

Semi Demountable Torch (HF)

[30-808-8039](#)

Suitable for ICPS-7500/8100 and ICPE-9000

Comprises:

Semi Demountable Torch Body [31-808-1258](#)

Semi Demountable Torch
 Adaptor [31-808-1260](#)

Tapered Alumina Injector 1.5mm [31-808-0013](#)

Joint Clip B14-Nickel [KCM-14](#)

GazFit Connectors (PKT 4)
 for 4mm OD side arms [GAZ-04](#)



HF Resistant Injectors

Tapered Alumina Injector 1.5mm [31-808-0013](#)

Tapered Alumina Injector 1.8mm [31-808-0015](#)

