

Eluo Nebulizer Cleaning Tool

As easy as 1, 2, 3

1. Fill with solvent



2. Insert nebulizer



3. Clean nebulizer



Particle build-up in a nebulizer capillary and tip causes sample flow to be constricted, thereby reducing nebulizer efficiency and performance. Now, blocked nebulizers can be safely and easily restored to optimum performance with the Eluo.

The Eluo is designed to efficiently deliver a cleaning solution through the nebulizer capillary to dislodge particle build-up and thoroughly clean the nebulizer. One simple action does it all. No more messy procedures or shattered nebulizers in ultrasonic baths. Use the Eluo regularly to maintain nebulizer performance and prolong nebulizer life. Every lab should have an Eluo.

Order P/N [70-ELUO](#) for glass concentric nebulizers or P/N [70-ELUO-OPD](#) for OpalMist, PolyCon or DuraMist nebulizers.

For more information visit: www.geicp.com

Ordering Information

Part Number	Description
70-ELUO	Eluo Nebulizer Cleaning Tool (For glass concentric nebulizers)
70-ELUO-OPD	Eluo Nebulizer Cleaning Tool (For OpalMist, PolyCon or DuraMist nebulizers)



P/N 70-ELUO-OPD



GLASS EXPANSION
Quality By Design

Head Office

6 Central Boulevard,
Port Melbourne, Vic 3207,
Australia
(61) 3 9320 1111
enquiries@geicp.com

Americas

31 Jonathan Bourne Drive,
Unit 7 Pocasset, MA 02559,
USA
508 563 1800
geusa@geicp.com

Europe

Friedenbachstrasse 9,
35781 Weilburg,
Germany
+49 6471 3778517
gegmbh@geicp.com