

FINALLY, A WORRY FREE ICP ARGON HUMIDIFIER

The Capricorn from Glass Expansion

The Argon Humidifier is commonly used in ICP analyses involving samples with high dissolved solids concentration to prevent salt build-up inside the sample introduction system.

The Glass Expansion Capricorn™ combines ease of operation, high performance and safety, to provide uninterrupted and maintenance-free operation.

- The rugged and durable plastic casing enhances safety.
- A built-in pressure indicator serves as a diagnostic so you always know when the Capricorn is under pressure.
- An easy-to-use bypass toggle allows you to take the Capricorn off-line without disconnecting argon lines.
- Maximum and minimum fill marks ensure that you are always operating under optimum conditions.
- A rugged bubbling frit ensures maximum saturation of the argon.
- The Capricorn is compatible with all ICP-OES and ICP-MS models.



INTERNATIONAL

Glass Expansion • 6 Central Boulevard • Port Melbourne • Vic 3207, Australia
Telephone: +61 3 9320 1111 • Toll Free (Aust): 1800 777 638
Facsimile: +61 3 9320 1112 • Email: enquiries@geicp.com

AMERICAS

Glass Expansion • 4 Barlows Landing Road • Unit 2A • Pocasset, MA 02559, USA
Telephone: 508 563 1800 • Toll Free (US): 800 208 0097
Facsimile: 508 563 1802 • Email: geusa@geicp.com

www.geicp.com



GLASS EXPANSION
Quality By Design